



LABOAO

Professional Laboratory Solution Provider

ZHENGZHOU LABOAO INSTRUMENT EQUIPMENT CO.,LTD



Web: www.laboao.com



Yd-1 Tablet Hardness Tester Manually

Model: YD-1





Tablet hardness testing, is a laboratory technique used by the pharmaceutical industry to test the breaking point and structural integrity of a tablet under conditions of storage, transportation, and handling before usage.

Video: https://www.youtube.com/embed/ojdxQYxrpRw

DESCRIPTION

Tablet hardness testers are instruments for detecting breaking hardness of tablet.



Email: sales@laboao.com Telephone: +86 371 61653992

Web: www.laboao.com

Applicable standards:

Industry standard (Tablet hardness tester)

FEATURES


- 1. The tablet is loaded and pressed artificially.
- 2. It has high accuracy press sensor and hardness data is displayed on LED.
- 3. It tests tablet hardness continuously.
- 4. Test data is displayed and locked automatically, so can test be reset and cycled.
- 5. Automation: auto-test, auto-diagnose, auto-alarm.

SPECIFICATIONS

Hardness range	2~296N 0.2kg~30kg
Hardness accuracy	±0.5%
Moving distance of probe	2~30mm
Power	220V/ 50Hz/ 10W or 110V/ 60Hz/ 10W
Dimension	320*200*100mm

Learn more about this product:

https://www.laboao.com/products/tablet-hardness-tester/yd-1-tablet-hardness-tester-manually



Choose LABOAO

Customize Your Exclusive Laboratory

High Quality

Perfect Service

Better Price

Find More











