

Fapas<sup>®</sup> website.

Fapas <sup>®</sup> QC MATERIAL DATA SHEET	T27327AQC
Matrix	Chocolate
Weight / Volume of Contents	20g

Analyte	Qualitative Results
Hazelnut or Hazelnut protein	NOT DETECTED (96% consensus)

This data sheet is applicable until	18 Jun 2025	
Recommended Storage on receipt Ambie		
Notes		
• Mix the QC material thoroughly before taking a representative analytical sample.		
<ul> <li>Stability of the QC material has been established as sufficient for the scope of the proficiency test from previous experience, expert advice and published literature. Fapas<sup>®</sup> advises that the QC material is analysed within the recommended date. Fapas<sup>®</sup> QC materials are intended to be used as single-analysis samples.</li> </ul>		
<ul> <li>Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the</li> </ul>		