



Fapas® QC MATERIAL DATA SHEET	T27335BQC
Matrix	Infant Formula (hydrolysed)
Weight / Volume of Contents	20g

Analyte	Qualitative Results
Milk	NOT DETECTED (93% consensus)

This data sheet is applicable until	11 Sep 2027
Recommended Storage on receipt	Ambient

Notes
<ul style="list-style-type: none">• Mix the QC material thoroughly before taking a representative analytical sample.• Stability of the QC material has been established as sufficient for the scope of the proficiency test from previous experience, expert advice and published literature. Fapas® advises that the QC material is analysed within the recommended date. Fapas® QC materials are intended to be used as single-analysis samples.• Full details on the proficiency test procedure used to characterise this QC material are available in the Protocol, Part 1 - Common Principles, freely available to download from the Fapas® website.